



## High Performance Covid Test SARS-COV-2 Antigen Rapid Test (For Covid 19) – Individually Wrapped

**Variant code: 22-COV19T/IW**

This listing is for 1 x box of Healgen Covid-19 Antigen Rapid Lateral Flow Tests. Each box contains 20 tests. The Healgen COVID-19 antigen rapid lateral flow test identifies infected individuals in just 15 minutes to increase workplace protection against the threat of COVID-19 and minimise disruption to businesses. It is ideal for working environments where staff are unable to carry out their job from home such as factories, retail premises, delivery services, healthcare service providers and hotels. Approved by the UK government and Public Health England, these rapid lateral flow tests are designed to detect large numbers of SARS-CoV-2 infectious cases in a short space of time to prevent further onward transmission. The rapid test is performed using nasopharyngeal swabs and is proven to give excellent accuracy (96.3%) with a sensitivity of 86.7% and a specificity of 100%. The test is simple to administer and comes supplied with clear instructions and all the necessary equipment to carry out the test. This includes a workstation to lay out the contents of the pack in a way that helps workplaces deliver mass testing with ease and efficiency to help keep them safe and operating effectively. Features Lateral flow test for rapidly detecting active infections Visual, easy to read result in just 15 minutes Tests can be administered anywhere, any time Helps to minimise disruption to businesses Ideal for businesses to regularly and reliably test staff Approved by the UK government and Public Health England In high demand from schools, councils and companies 96.3% accuracy High sensitivity and specificity Sold in a pack of 20 for testing en masse All necessary reagents provided and no additional equipment needed Each kit includes clear instructions and is easy to use 24 month shelf life Storage between 2 and 30°C CE certified It is understood that as many as one in three people who have coronavirus do not show symptoms but are infectious. Lateral flow tests help to quickly identify and reduce the spread of the virus from those who appear healthy. For further details, EC Declaration of Conformity and usage instructions, please see 'Product

Downloads'. Pack Contents: 20 Test Cassettes 20  
Sterile Swabs 20 Extraction Tubes and dropper  
tips 1 Workstation 2 Buffers 1 Package Insert