



Evans Protect Disinfectant Cleaner Concentrate Refill – 2 x 5L

Variant code: 05-215

Evans Protect is a professional **disinfectant cleaner concentrate** developed by Evans Vanodine for use in hygiene-critical environments. This multipurpose formulation is designed to clean and disinfect in one application, helping maintain high standards of cleanliness and infection control across a wide range of surfaces.

Suitable for use in environments such as **schools, hospitals, surgeries, nursing homes, sports clubs and leisure facilities**, Protect delivers effective bactericidal performance while removing general soiling. It is particularly suited to the leisure industry for cleaning equipment, sunbeds, steam rooms and saunas. The formulation passes **BS EN 1276** standards and is effective against a wide range of bacteria, including **MRSA and Shigella sonnei**, as well as viruses and fungi, helping to reduce the risk of cross infection. It leaves a **clean, fresh odour** after use and is suitable for a variety of hard surfaces.

Product Features

- **Multipurpose disinfectant cleaner concentrate**
- Cleans and disinfects in **one application**
- Effective against a **wide range of bacteria, viruses and fungi**
- Passes **BS EN 1276** hygiene standard
- Suitable for a variety of **hard surfaces**
- Ideal for use in **healthcare and leisure environments**
- Leaves a **clean, fresh fragrance**

Product Benefits

- Supports infection control in high-risk environments
- Reduces risk of cross-contamination on frequently touched surfaces
- Improves efficiency by combining cleaning and disinfection
- Provides versatile use across multiple sectors
- Maintains hygiene standards in demanding environments

- Enhances cleanliness and freshness of cleaned areas
- Suitable for routine and deep cleaning applications

Technical Characteristics

- Designed for use as a **concentrate requiring dilution**
- Suitable for use on a wide range of **hard washable surfaces**
- Effective against bacteria, viruses and fungi when used as directed
- Passes recognised **British and European test standards**
- Apply using appropriate cleaning methods depending on surface and environment
- Follow manufacturer instructions for dilution and contact time
- Store in a cool, dry environment away from direct sunlight

Applications

Evans Protect is suitable for use in healthcare facilities, educational environments, leisure centres, sports clubs and general commercial settings where infection control is essential. It is ideal for cleaning and disinfecting surfaces such as equipment, work areas and high-touch points, helping maintain hygienic conditions and reduce the spread of harmful microorganisms in professional environments.

Property	Value
Fragrance	Floral
Usage Instructions	Suitable for use on all hard wood washable surfaces
Dilution	Dilute with clean water – viricidal application 50ml per 750ml diluted solution. Bactericidal applications 30ml per 750ml diluted solution
pH	12
CLP	corrosive, irritant
How to Use	Clean surfaces before disinfection. apply to surface by spray, clean cloth or mop. Leave on surface for required contact time. Wipe with cloth and air dry
Chemical Colour	Blue
EN1276 Certified	Yes
EN1650 Certified	Yes

Property	Value
General Cleaning	For non-disinfection cleaning tasks. Dilute 1:50 with clean water (100 ml per 5 L diluted solution). Mop or wipe surface.
EN16777 Certified	Yes
EN16615 Certified	Yes
Showerhead Cleaning	Soak showerhead in a solution of Protect at a dilution rate of 1:100 for a minimum of 5 minutes. This procedure should be carried out every 3 months.