



## Milton Disinfecting Liquid – 5L

**Variant code: 04-250**

Milton Disinfectant is a powerful and safe disinfecting fluid with a variety of professional uses. The chlorine-based, concentrated formula effectively kills bacteria, fungi and viruses in under 15 minutes, making it ideal for use in hospitals and nursing homes.

Milton leaves no trace of taste or smell, so it is also suitable to be used within food preparation and is also safe enough to be used to sterilise baby bottles, teats and feeding equipment.

1 litre of concentrated fluid makes 22 litres of solution. This 5L bottle provides enough disinfecting fluid to make more than 200 mop buckets or 450 spray bottles.

Conforms to European Standards EN1276 and EN13697.

### Uses Include

- Diluting into mop bucket or spray bottle for sanitising hard surfaces
- Sterilising baby bottles and teats
- Cleaning bed pans and tools in healthcare environments
- Sanitising water before drinking
- Washing food and utensils in food preparation environments

### How to Use Milton Disinfectant

- Ensure surface dirt is cleaned off before sterilising
- Prepare mop bucket or spray bottle with appropriate amount of water
- Add correct dose of Milton Disinfecting Fluid as per instructions on packaging
- Submerge items or spray onto surface to be disinfected
- Leave for 15 minutes
- Dry items ready for use – no need to rinse!
- Items and surfaces will remain sterilised for 24 hours

**Property**

**Value**

Property	Value
Fragrance	Characteristic
EN	EN1276, EN13697
Dilution	As per instructions on packaging
pH	11.8
Area of Use	Floors, healthcare, drinking water, food preparation areas and equipment
CLP	irritant, environmental
How to Use	Remove soiling prior to sterilising Prepare mop bucket or spray bottle. Submerge items or spray onto surface to be disinfected. Leave for 15 minutes. Dry items ready for use, no need to rinse.
Chemical Colour	Yellow